


General Combustion Supplies

LECO#	Sci-Lab#	Description	UOM
774-649	AR653	Wire Cleaning Brush Assembly	EA
774-204	AR4204	Stainless Steel Crucibles	10
502-461	AR2461	Cotton Thread Fuse	375
502-462	AR205	Fuse Wire	375
502-460	AR02460	Fixed Wire Ignitor	5



Reference Materials


LECO#	Sci-Lab#	Description	UOM
502-892	AR208C	Benzoic Acid Tablets, 1g	100
	AR208	Benzoic Acid Tablets, 1g	10
	AR208V	Benzoic Acid Tablets, 1g	500


O-rings

LECO#	Sci-Lab#	Description	UOM
611-543	AR1543	O-ring	EA
621-605-560	AR5560	O-ring	EA
621-605-664	AR5664	O-ring	EA
772-738	AR738	O-ring	EA
773-913	AR3913	O-ring	EA
621-656-001	AR6001	O-ring	EA
608-624	AR624	O-ring - AC500	EA
608-750	AR750	O-ring - AC500	EA
760-224	AR224	O-ring - AC500	EA

Operational Items

LECO#	Sci-Lab#	Description	UOM
501-241	AR241	Silicon Vacuum Grease	35g
621-657	AR21657	Filter Cartridge, 1 Micron	EA
776-978	AR978	Large Capacity Glass Scoop	EA
502-051	AR1618	Tweezers, Stainless Steel, Curved Point, 4.5" Length	EA





Met Coke
AR2771



Coal
AR2773

The Ultimates (50 grams)

Sci-Lab#	%Ash	%Vol.	%Fixed C	BTU	%Sulfur	%Carbon	%Hydrogen	%Nitrogen	%Chlorine	%Oxygen
AR2771 Met Coke	7.53	(0.79)	(91.68)	13178	0.57	89.89	(0.23)	1.08	(0.024)	(0.70)
AR2772 Met Coke	9.19	(1.0)	(89.82)	12777	0.80	88.39	(0.15)	1.14	(0.029)	(0.33)
AR2773 Coal	5.03	29.48	65.49	12572	0.52	79.7	3.13	0.64	(<0.01)	10.98
AR2775 Coal	6.31	43.54	(50.15)	11888	0.41	69.69	5.09	1.00		(17.5)
AR2776 Coal	21.02	15.76	(63.22)	11232	0.81	68.88	2.76	1.03	(0.0385)	(5.5)
AR2778 Coal	5.54	17.58	(76.88)	14829	0.76	86.20	4.36	1.20		(1.94)
AR2780 Coal	19.33	30.37	(50.30)	11571	3.04	65.03	4.38	1.23	(0.144)	(6.99)
AR2781 Coal	42.89	26.30	(30.81)	8097	2.27	45.14	3.39	0.91		(5.40)
AR2782 Coal	23.64	32.24	44.12	10031	5.39	55.70	4.10	(1.28)	(0.0282)	(9.89)
AR2783 Anthracite	20.97	7.01	(72.02)	11336	0.54	72.32	2.03	0.90	(0.0224)	(3.24)

Sci-Lab Analytical Ltd. Is a the exclusive UK dealer for Alpha Resources. Alpha Resources provides industry leading references materials with fit-for-purpose results. We work with Alpha Resources for the assurance of quality that comes with an ISO 9001, 17025 and 17034 accredited supplier.

Did you know Alpha's lab certification process was featured in Lab Manager (April 2021)? [Click here for details](#)

Alpha Resources offers an extensive array of oil CRMs. For a comprehensive list of Alpha's in-stock and available oil CRMs, click here to see their website for a full list. We have an extensive warehouse of Alpha products in the UK, and receive shipments weekly. You can be assured of fast delivery of your lab supplies.