


Quartz, Glass, and Furnace Supplies

Elementar#	Sci-Lab#	Description	UOM
11.00-1170/4	AR1170/4	Ash Finger	EA
11.00-1170/4	AEC3060	Ash Finger - 10 Pack	EA
	AEC1020	Cap Kit For Threaded Joints Capts PTFE Backed Silicone Seals	2
11.00-1150/4	AR11504	Combustion/Reduction Tube	EA
22 160 200/4	AEC3008	Drying Tube	EA
15.00-0080	AEC1143	Drying Tube	EA
22 160 204/4	AEC3066	Drying Tube, Straight	EA
11001317/4	AR1317/4	Protection tube Vario (105mm) (with slits)	EA
11001317/4	AEC3059	Protection tube Vario (105mm) (without slits)	EA
15.40-1002	AEC3057	Support Rod Reduction Tube (110mm)	EA
23.00-1148	AEC3052	Support Rod Reduction Tube (175mm)	EA
11.00-1364/4	AEC3053	Support Rod Reduction Tube (65mm)	EA
11.00-1172/4	AR1172/4	Quartz Support Tube	EA
11 001 171/4	AR11714	Support Tube (133mm)	EA
	AEC3123	U-Tube	EA





Reference Card: *Elementar*[®] UNICUBE/vario MICRO cube



Operational Accessories


Elementar#	Sci-Lab#	Description	UOM
200007948/ 11.00-1363/4	AR1363/4	Ash Crucible, Steel	EA
200004147/ 11.00-1331/4	AR13314	Connecting Tube Transparent Silica Includes O-rings	EA
	AEC3006	Oxygen Lance (130mm)	EA
	AEC3006- 10	Oxygen Lance (130mm)	10
200014764/ 05 000 088	AR50088	Oxygen Lance (140mm)	EA
	AR50088- 10	Oxygen Lance (140mm)	10
200004159/ 05 000 354	AR00354	Oxygen Lance (150mm)	EA
200010997/ 41.01-0003/4	AEE2535	Manual Press (7mm Pellet)	EA
200011004/ 15.00-0005	AEE2532	Set of Ball Valve Half Shells (10mm)	EA
200006998/ 03 968 006	AR241	Silicone Grease	35g
200013169/ 05 001 715	AR01715	Cotton Wool	50g
200015256/ 03 002 399	AEE1174	Filter (18mm)	1
	AEE1174-10	Filter (18mm)	10
200006785/ 05 000 094	AEE2505	Gloves	Pair
200009657/ 11.00-0500/4	AEE2521	Manually Operated Drill	EA







Chemical Reagents

Elementar#	Sci-Lab#	Description	UOM
	AR01008-200	Alumina (Aluminium Oxide) Spheres (3-6mm)	200g
05 008 467	AR01008	Alumina (Aluminium Oxide) Spheres (3-6mm)	40g
05 000 086	AR00086	Alumina (Aluminium Oxide) Wool Fine	20g
05 000 699	AR189	Copper Sticks	100g
	AEB1206	Copper Fine Wires Reduced (4 x 0.5mm)	1kg
05 001 039	AR01039	Copper Oxide Wire Fine Wires (4 x 0.5mm)	250g
	AEB1209	Copper Oxide Wire Fine Wires (4 x 0.5mm)	3kg
	AEB1091M	Copper Wires Fine Wires Reduced (4 x 0.5mm)	3 x 100g
	AEB1083M	Copper Wire Fine Wires Reduced (4 x 0.5mm)	454g
03 679 900	AR168	AlphaSolve® II	250g
	AEB1255	Quartz Wool Discs (24mm)	25
03 679 907	AEB1008	Silica (Quartz) Chips Granular (1-2mm)	100g
03 679 908	AR177	Quartz Wool, Fine	50g
	AEB1276	Silica (Quartz) Wool Very Fine	15g
22 131 365	AR22365	Silver Wool, Fine Wire (0.05mm)	50g
11.02-0008	AR2181	Tungsten (VI) Oxide Granular	60g
11.02-0024	AR20024	Tungsten (VI) Oxide Powder	90g
11.02-0017	AR20017	Tungsten (VI) Oxide Powder Fine Powder	25g




Prepacked Reactors

Elementar#	Sci-Lab#	Description	UOM
	AEX2505	Prepacked CNS/CHNS Reaction Tube for Micro Cube	EA
	AEX2506	Prepacked N/NC/CHN Reaction Tube Vario for EL3/cube/Micro Cube	EA
	AEX3026	Prepacked Filled Trap	EA
	AEX3022	Prepacked Filled Trap	EA



O-rings

Elementar#	Sci-Lab#	Description	UOM
200014191/ 05 001 302	AEE1087	Kalrez O-ring Black 7.59 x 2.62mm	EA
	AEE1021	O-ring Schott Type Silicone PTFE	10
200015397/ 05 000 369	AR00369	O-ring (4 x 4mm)	EA
200013028/ 05 000 096	AR00096	O-ring (18 x 2mm)	EA
200007224/ 03 654 627	AR54627	O-ring (20 x 2mm)	EA
	AR54627-10	O-ring (20 x 2mm)	10
200007469/ 03 002 267	AR02267	O-ring (24 x 2mm)	EA
200007469/ 03 002 267	AR02267-10	O-ring (24 x 2mm)	10
200008314/ 05 000 249	AR00249	O-ring (4.5 x 1.5mm)	EA
200008314/ 05 000 249	AR00249-10	O-ring (4.5 x 1.5mm)	10
200006004/ 05 000 416	AR00416	O-ring (7.5 x 2.5mm)	EA
200006004/ 05 000 416	AR00416-10	O-ring (7.5 x 2.5mm)	10
200005999/ 05 000 353	AR50353	Quadring (7.65 x 1.78mm)	EA
200005999/ 05 000 353	AR50353-10	Quadring (7.65 x 1.78mm)	10
200010515/ 05 000 247	AR00247	Seal for Oxygen Lance Rubber	EA



Aluminium Boats

Elementar#	Sci-Lab#	Description	UOM
03 674 012	AED5132	Aluminium Weighing Boat (12 x 4 x 4mm)	1000
03 674 016	AED5133	Aluminium Weighing Boat (12 x 6 x 6mm)	500

Aluminium Capsules for Liquids

Elementar#	Sci-Lab#	Description	UOM
05 000 814	AR00814	Aluminium Capsules, Smooth Wall, Flat Base (7 x 2.7mm)	250

Silver Boats

Elementar#	Sci-Lab#	Description	UOM
	AED5127	Silver Weighing Boat (6 x 4 x 4mm)	500
	AED5024	Silver Weighing Boat (12 x 4 x 4mm)	100
11.00-1371	AR01371	Silver Weighing Boat (12 x 4 x 4mm)	500
22 133 212	AR33212	Silver Weighing Boat (12 x 6 x 6mm)	100
22 133 213	AR33213	Silver Weighing Boat (15 x 8 x 8mm)	70
22 133 214	AR33214	Silver Weighing Boat (20 x 10 x 9mm)	42

Silver Capsules for Liquids

Elementar#	Sci-Lab#	Description	UOM
05 000 816	AR00816	Silver Capsule, Smooth Wall (5.5 x 3mm)	250
05 001 732	AED4103	Silver Capsule, Smooth Wall (5.5 x 2mm)	250
	AED4101	Silver Capsule, Smooth Wall (9 x 3.5mm)	100
05 000 817	AED4109	Silver Capsule, Smooth Wall (9 x 3.5mm)	250

Tin Boats

Elementar#	Sci-Lab#	Description	UOM
22 137 418	AR37418	Tin Boat (12 x 4 x 4mm)	1000
22 137 419	AR37419	Tin Boat (12 x 6 x 6mm)	500
22 137 420	AR37420	Tin Boat (15 x 8 x 8mm)	250
12.00-1009	AR01009	Tin Boat (20 x 10 x 9mm)	120
11.00-5057	AED5140	Tin Boat Heavy Weight (12 x 4 x 4mm)	1000
11.02-1047	AED5029	Tin Weighing Boat (6 x 4 x 4mm)	1000

Tin Cones Squares and Discs

Elementar#	Sci-Lab#	Description	UOM
05 000 429	AR2186	Tin foil Cones, Standard Weight	100
05 002 797	AR2186-1000	Tin foil Cones, Standard Weight	1000
12.01-0032	AR10032	Tin Foil Squares, Standard Weight (37 x 37mm)	400
12.01-0033	AR2397-400	Tin Foil Squares, Standard Weight (50 x 50mm)	400

Tin Capsules for Liquids

Elementar#	Sci-Lab#	Description	UOM
05 000 481	AR0481	Tin Capsule, Smooth Wall, Flat Base (11 x 4.75mm)	250
05 000 312	AED4120	Tin Capsule, Smooth Wall, Flat Base (14 x 6.75mm)	100
05 001 727	AED4033	Tin Capsule, Smooth Wall, Flat Base (5.5 x 3.5mm)	100
	ATD4001	Tin Capsule, Smooth Wall, Flat Base (6 x 3mm)	100
03 951 620	AED4131	Tin Capsule, Smooth Wall, Flat Base (9 x 3mm)	250
	AED4035	Tin Capsule, Smooth Wall, Flat Base (9 x 3.5mm)	100
03 951 619	AED4132	Tin Capsule, Smooth Wall, Flat Base (9 x 3.5mm)	250
05 000 311	AR2166	Tin Capsule, Smooth Wall, Round Base (20 x 8mm)	100

Tin Capsules/Pans for Solids

Elementar#	Sci-Lab#	Description	UOM
05 003 400	AED1103	Tin Capsule, Pressed, Standard Weight (11 x 8mm)	100
	ATD1010	Tin Capsule, Pressed, Standard Weight (12.5 x 5mm)	200
05 003 395	ATD1002	Tin Capsule, Pressed, Standard Weight (5 x 3.5mm)	500
	ATD1009	Tin Capsule, Pressed, Standard Weight (8 x 5mm)	100
05 003 394	ATD1008	Tin Capsule, Pressed, Standard Weight (8 x 5mm)	250
	ATD1021	Tin Capsule, Pressed, Standard Weight (9 x 5mm)	250
	ATD1022	Tin Capsule, Pressed, Standard Weight (9 x 5mm)	400
	ATD1013	Tin Capsule, Pressed, Ultra-Light Weight (8 x 5mm)	100
	AED1019	Tin Containers, Light Weight (8 x 5mm)	1000