

## Reagents, Catalyst, Oxidation, and Supplies

Elementar#	Sci-Lab#	Description	UOM
	AEB1083M	Copper Wires, Fine Wires Reduced, Pound Pack (4 x 0.5mm)	454g
200010245/ 50 008 467	AEB1090	Alumina (Aluminium Oxide) Spheres (3-6mm)	50g
	AEB1093	Alumina (Aluminium Oxide) Spheres (3-6mm)	250g
200007811/ 03 002 262	AEB1196	Platinum Catalyst, 5% Pt	30g
	AEB1206	Copper Oxide Wires, Fine Wires (4 x 0.5mm)	1kg
	AEB1209	Copper Oxide Wires, Fine Wires (4 x 0.5mm)	3kg
	AEB1246	Platinum Catalyst, 0.5% Pt	30g
200010160/ 38.00-0124	AEB1304	Brass Wool	50g
200008093/ 16 02 0056	AEB1305	Post Combustion Reagent for Exceed® Instruments	500g
200007811/ 03 002 262	AEB1308	Platinum Catalyst, 5% Pt	150g
200004773/ 05 000 086	AR00086	Alumina Wool, Fine	20g
200009554/ 05 001 039	AR01039	Copper Oxide Wires, Fine Wires	250g
200004827/ 03 679 908	AR177	Quartz Wool, Fine	50g
200014271/ 05 000 699	AR189	Copper Sticks	100g
200005872/ 22 131 365	AR22365	Silver Wool, Fine Wire (.05mm)	50g



## Reference Card: *Elementar®* *rapid N exceed*



## Tin Capsules for Liquids

Elementar#	Sci-Lab#	Description	UOM
200013285/ 05 000 312	AED4120	Tin Capsules, Smooth Wall (14 x 6.75mm)	100
200008184/ 05 000 481	AR0481	Tin Capsules, Pressed (11.5 x 4.75mm)	250


## Tin Cones, Squares, and Discs

Elementar#	Sci-Lab#	Description	UOM
200009010/ 05 000 429	AR2186	Tin Foil Cups (36 x 36mm)	100
200009795/ 05 002 797	AR2186-1000	Tin Foil Cups (36 x 36mm)	1000
	AR2397	Tin Foil Cups (51 x 51mm)	100
	AR2397-400	Tin Foil Cups (51 x 51mm)	400

## Quartz, Glassware, and Furnace Tubes

Elementar#	Sci-Lab#	Description	UOM
200011489/ 16.00-1012	AEC3026	Drying Tube (348mm)	EA
200010584/ 16.00-1037	AEC3047	Drying Tube 150mm (with GL25 Caps)	EA
	AEC3049	Quartz Reduction/Tungsten Tube modified with dimple (29mm)	EA
200006484/ 16.02-0036	AEC3069	Combustion Tube	EA
200013679/ 25.00-1048	AR01048	Post Combustion Tube	EA





## Supplies for Furnace Tubes

Elementar#	Sci-Lab#	Description	UOM
	AEC3006	Oxygen Lance (130mm)	EA
	AEC3006-10	Oxygen Lance (130mm)	10
200004619/ 16.00-1104	AEC3065	Ash Crucible, Stainless Steel	EA
200011310/ 12.00-1004/4	AR01004/ 1-10	Wire Mesh Ronde (22mm)	10
200011310/ 12.00-1004/4	AR01004/4	Wire Mesh Ronde (22mm)	EA
200010457/ 12.01-1008/4	AR02530S	Tungsten Spacers, for use with AR02530	25
	AR02530S-100	Tungsten Spacers, for use with AR02530	100
200011961/ 12.01-1032/4	AR11032/4	Distance Ring	EA
200004147/ 11.00-1331/4	AR13314	Connecting Tube Transparent Silica (65mm) (includes O-rings)	EA
200014764/ 05 000 088	AR50088	Oxygen Lance (140mm)	EA
200014764/ 05 000 088	AR50088-10	Oxygen Lance (140mm)	10

## Other Sample Encapsulation

Elementar#	Sci-Lab#	Description	UOM
200006862/ 16.02-0055	AEX2507	Prepacked reduction tube Exceed® Instruments, Quartz	EA

## Operational Accessories

Elementar#	Sci-Lab#	Description	UOM
	AEE2001	Scoop, SS Ladle (300mm)	EA
200006785/ 05 000 094	AEE2505	Gloves, Temperature Resistant	Pair
200009514/ 05 001 432	AEE2520	Brush Brush	EA
200009657/ 11.00-0500/4	AEE2521	Manually Operated Drill for Cleaning Tubes	EA
200004971/ 41.01-0004	AEE2536	Manual Pellet Press (9mm)	EA
200009549/ 41.01-0002/4	AEE2537	Manual Pellet Press (10.5mm)	EA
200013169/ 05 001 715	AR01715	Cotton Wool	50g
200010676/ 03 001 248	AR1618	Tweezers, Stainless Steel, Curved Point (4.5in)	EA
200014862/ 03 968 006	AR241	Silicone Grease	35g

## Prepacked Tubes

Elementar#	Sci-Lab#	Description	UOM
200006862/ 16.02-0055	AEX2507	Prepacked reduction tube Exceed Instruments Quartz	EA

## Prepacked Tubes

Elementar#	Sci-Lab#	Description	UOM
200007904/ 05 000 568	AEE1175	O-ring (16 x 3mm)	EA
200007904/ 05 000 568	AEE1175-10	O-ring (16 x 3mm)	10
200015449/ 05 002 707	AEE1177	O-ring (25 x 2mm)	EA
200015449/ 05 002 707	AEE1177-10	O-ring (25 x 2mm)	10
200007217/ 05 000 095	AR00095	O-ring (22 x 2mm)	EA
200007217/ 05 000 095	AR00095-10	O-ring (22 x 2mm)	10
200013028/ 05 000 096	AR00096	O-ring (18 x 2mm)	EA
200010515/ 05 000 247	AR00247	Seal for Oxygen Lance - 10cm piece which can be cut into ~30 x 3mm seals	10
200008314/ 05 000 249	AR00249	O-ring (4.5 x 1.5mm)	EA
200008314/ 05 000 249	AR00249-10	O-ring (4.5 x 1.5mm)	10
200008183/ 05 000 0359	AR00359	O-ring (9.19 x 2.62mm)	EA
200008183/ 05 000 0359	AR00359-10	O-ring (9.19 x 2.62mm)	10
200006004/ 05 000 416	AR00416	O-ring (7.5 x 2.5mm)	EA
200006004/ 05 000 416	AR00416-10	O-ring (7.5 x 2.5mm)	10
200005999/ 05 000 353	AR50353	Square Gasket (7.65 x 1.78mm)	EA
200005999/ 05 000 353	AR50353-10	Square Gasket (7.65 x 1.78mm)	10
200007224/ 03 654 627	AR54627	O-ring (20 x 2mm)	EA
200007224/ 03 654 627	AR54627-10	O-ring (20 x 2mm)	10