



Quartz, Glassware, and Furnace Tubes


Elementar#	Sci-Lab#	Description	UOM
200005357/ 25.00-1028	AR01028	Combustion Tube, Quartz, vario MAX	EA
200014722/ 25.00-1030	AR01030	Combustion Tube vario MAX, Stainless Steel	EA
200013540/ 22 160 200/4	AEC3008	Drying Tube	EA
200013540/ 22 160 204/4	AEC3066	Drying Tube, Straight	EA
200013679/ 25.00-1048	AR01048	Post Combustion Tube	EA
200008838/ 25.00-0017	AR01007	Quartz Combustion Tube	EA
200013606/ 34.00-1339/4	AEC3038	U-Tube Large Complete with GL18 Caps	EA
200006654/ 25.00-1123	AEC3058	Ceramic Crucibles vario MAX	100
00014515/ 03 002 278	AR2278	Metal Crucibles	100



Operation Accessories

Elementar#	Sci-Lab#	Description	UOM
200004147/ 11.00-1331/4	AR13314	Connecting Tube Transparent Silica includes O-rings	EA
200011961/ 12.01-1032/4	AR11032/4	Distance Ring	EA
200005508/ 25.00-0029	AR00029	Elementar Gripper Arm Assembly, with AR02395 O-ring	EA
200010457/ 12.01-1008/4	AR02530S	Tungsten Spacers for Elementar, for use with AR02530	25
200011310/ 12.00-1004/4	AR01004/4	Wire Mesh Ronde (22mm)	EA
200006395/ 50 008 449/4	AR08449/4	Wire Screen (26.8mm)	EA
200009514/ 05 001 432	AEE2520	Brass Brush	EA
200013169/ 05 001 715	AR01715	Cotton Wool	50g
200015256/ 03 002 399	AEE1174	Filter Diameter (18mm)	1
200015256/ 03 002 399	AEE1174-10	Filter Diameter (18mm)	10
200014862/ 03 968 006	AR241	Silicone Grease	35g





Chemical Reagents

Elementar#	Sci-Lab#	Description	UOM
	AEB1093	Alumina (Aluminium Oxide) Spheres (3-6mm)	250g
200010245/ 50 008 467	AR01008	Corundum Balls	40g
200014271/ 05 000 699	AR189	Copper Sticks	100g
	AR33214	Silver Sample Boats (20 x 10 x 9mm)	42
	AEB1206	Copper Oxide Wires Fine Wires (4 x 0.5mm)	1kg
200009554/ 05 001 039	AR01039	Copper Oxide Wire, Fine Wires	250kg
	AEB1209	Copper Oxide Wires, Fine Wires (4 x 0.5mm)	3kg
	AEB1091	Copper Wires Pack, Fine Wires Reduced (4 x 0.5mm)	3 x 100g
	AEB1083M	Copper Wires, Fine Wires Reduced, Pound Pack (4 x 0.5mm)	454g
200013014/ 05 005 080	AR171	Magnesium Perchlorate	454g
200011868/ 03 679 900	AR168	AlphaSolve® II	250g
	AEB1246	Platinum Catalyst 0.5% Pt	30g
200007811/ 03 002 262	AEB1196	Platinum Catalyst 5% Pt	30g
	AEB1255	Quartz Wool discs (24mm)	25
200007521/ 05 000 900	AR00900	Quartz Chips (Turnings), (2-5mm)	50g
200004827/ 03 679 908	AR177	Quartz Wool, Fine	10g
	AEB1276	Silica (Quartz) Wool, Very Fine	15g
200005872/ 22 131 365	AR22365	Silver Wool, Fine Wire (.05mm)	50g
200010974/ 12.00-0040	AR02530	Tungsten Metal	500g
200010174/ 12.01-0045	AR02530X	Tungsten Kit Complete with 500G (AR02530) of Tungsten	KIT



Corundum Balls
AR01008



Copper Sticks
AR189

O-ring

Elementar#	Sci-Lab#	Description	UOM
200010589/ 03 002 395	AR02395	O-ring for Gripper Arm	EA
03 654 701	AR54701	O-ring Silicone (4 x 4mm)	EA
200007224/ 03 654 627	AR54627	O-ring (20 x 2mm)	EA
200007224/ 03 654 627	AR54627-10	O-ring (20 x 2mm)	10
200007217/ 05 000 095	AR00095	O-ring (22 x 2mm)	EA
200007217/ 05 000 095	AR00095-10	O-ring (22 x 2mm)	10
200008314/ 05 000 249	AR00249	O-ring (4.5 x 1.5mm)	EA
200008314/ 05 000 249	AR00249-10	O-ring (4.5 x 1.5mm)	10
200006004/ 05 000 416	AR00416	O-ring (7.5 x 2.5mm)	EA
200006004/ 05 000 416	AR00416-10	O-ring (7.5 x 2.5mm)	10
200008183/ 05 000 359	AR00359	O-ring (9.19 x 2.62 mm)	EA
200008183/ 05 000 359	AR00359-10	O-ring (9.19 x 2.62 mm)	10
200008068/ 03 002 268	AR02268	O-ring (32 x 3mm)	EA
200006094/ 05 000 425	AR00425	O-ring (5 x 1.5mm)	EA
200009565/ 03 002 850	AEE1176	O-ring (4 x 2mm)	EA
200005999/ 05 000 353	AR50353	Quading (7.65 x 1.78mm)	EA
200005999/ 05 000 353	AR50353-10	Quading (7.65 x 1.78mm)	10
200006143/ 03 002 607	AR02607	Shaft Seal	EA